



---

## PRODUCT INFORMATION

---

### **Benchtop External Cavity Laser Controller**

#### **Laser Current | TEC | Piezo Scanning | Remote Control**

The PilotPZ is a laser diode controller together with laser diode temperature controller, Piezo Amplifier and Function Generator combined in one bench-top instrument. Special features are low noise laser current and highly stable temperature control. Remote control via USB, GPIB / IEEE 488.2, and RS232 interface.

---

## CONTACT & LEGAL

---



For further information about the product or support requests, please contact your local person in charge.

+ 49 (0)6421 305-0 Germany  
+ 1 800 352 3639 United States, Canada

Sacher Lasertechnik GmbH  
Rudolf-Breitscheid-Str. 1-5  
35037 Marburg  
Germany

Sacher Lasertechnik LLC.  
5765 Equador Way  
Buena Park, CA 90620  
United States



Copyright Sacher Lasertechnik Group 2010  
<http://www.sacher-laser.com>  
All Rights Reserved.

Sacher Lasertechnik and the Sacher Lasertechnik logo are trademarks or registered trademarks of Sacher Lasertechnik GmbH - Germany, other countries, or both.

The Lasers Systems Lynx™, Tiger™, Cheetah™, Lion™, Cougar™, Jaguar™, Serval™, Serval Plus™ and CAT™ are trademarks or registered trademarks of Sacher Lasertechnik GmbH - Germany, other countries, or both.

Other company, product and service names may be trademarks or service marks of others. References in this publication to Sacher Lasertechnik products or services do not imply that Sacher Lasertechnik intends to make them available in all countries in which Sacher Lasertechnik operates.

