



---

## PRODUCT INFORMATION

---

### **Motorized Tunable Littman/Metcalf Laser System - X500**

**397 nm .. 2450 nm | up to 100 mW | motor tuning | patented technology**

Motorized Littman / Metcalf Motorized Tunable External Cavity Diode Laser, 630nm .. 2450nm, up to 100mW, hands-off operation, fully remote controlled wavelength scans, large mode-hop free piezo tuning range, narrow linewidth, free space and fiber coupled versions. Check our ar-coated diode stock list for availability.

---

## CONTACT & LEGAL

---



For further information about the product or support requests, please contact your local person in charge.

+ 49 (0)6421 305-0 Germany  
+ 1 800 352 3639 United States, Canada

Sacher Lasertechnik GmbH  
Rudolf-Breitscheid-Str. 1-5  
35037 Marburg  
Germany

Sacher Lasertechnik LLC.  
5765 Equador Way  
Buena Park, CA 90620  
United States



Copyright Sacher Lasertechnik Group 2010  
<http://www.sacher-laser.com>  
All Rights Reserved.

Sacher Lasertechnik and the Sacher Lasertechnik logo are trademarks or registered trademarks of Sacher Lasertechnik GmbH - Germany, other countries, or both.

The Lasers Systems Lynx™, Tiger™, Cheetah™, Lion™, Cougar™, Jaguar™, Serval™, Serval Plus™ and CAT™ are trademarks or registered trademarks of Sacher Lasertechnik GmbH - Germany, other countries, or both.

Other company, product and service names may be trademarks or service marks of others. References in this publication to Sacher Lasertechnik products or services do not imply that Sacher Lasertechnik intends to make them available in all countries in which Sacher Lasertechnik operates.

