



PRODUCT INFORMATION

Tapered Littrow Laser - Tiger

750nm .. 1080nm | up to 1000mW manual tuning | patented technology

Monolithic Tapered Littrow Tunable External Cavity Laser, 750nm .. 1080nm, up to 1000mW, patented technology, monolithic cavity design, hands-off operation, narrow linewidth, up to 30GHz mode-hop free tuning, free space and fiber coupled versions. Check our tapered diode stock list for availability.

CONTACT & LEGAL



For further information about the product or support requests, please contact your local person in charge.

+ 49 (0)6421 305-0 Germany
+ 1 800 352 3639 United States, Canada

Sacher Lasertechnik GmbH
Rudolf-Breitscheid-Str. 1-5
35037 Marburg
Germany

Sacher Lasertechnik LLC.
5765 Equador Way
Buena Park, CA 90620
United States



Copyright Sacher Lasertechnik Group 2010
<http://www.sacher-laser.com>
All Rights Reserved.

Sacher Lasertechnik and the Sacher Lasertechnik logo are trademarks or registered trademarks of Sacher Lasertechnik GmbH - Germany, other countries, or both.

The Lasers Systems Lynx™, Tiger™, Cheetah™, Lion™, Cougar™, Jaguar™, Serval™, Serval Plus™ and CAT™ are trademarks or registered trademarks of Sacher Lasertechnik GmbH - Germany, other countries, or both.

Other company, product and service names may be trademarks or service marks of others. References in this publication to Sacher Lasertechnik products or services do not imply that Sacher Lasertechnik intends to make them available in all countries in which Sacher Lasertechnik operates.

